


## XRF Platinum Program for Herzog



### XRF Crucibles suitable for Herzog (without rim)




Diameter internal		Height (H)	Base thickness (s)	Weight approx.	Heraeus drawing
top (OD)	bottom (UD)				
mm	mm	mm	mm	g	
50	38,5	30	0,8	76	80047806
50	38,5	30	1,0	90	80047806



### XRF Crucibles suitable for Herzog (with rim)




Diameter internal		Height (H)	Base thickness (s)	Weight approx.	Heraeus drawing
top (OD)	bottom (UD)				
mm	mm	mm	mm	g	
50	38,5	30	1	115	80303231




### XRF Crucibles suitable for Herzog (with crucible holder)





Diameter internal		Height (H)	Base thickness (s)	Weight approx.	Heraeus drawing
top (OD)	bottom (UD)				
mm	mm	mm	mm	g	
36	22	38	0,6	64	80623511






#### XRF Casting dishes suitable for Herzog (square version)

Diameter internal		Flange Length (FL)	Height (H)	Base thickness (s)	Weight approx.	Heraeus drawing
top (OD)	bottom (UD)					
mm	mm	mm	mm	mm	g	
41	38	50	4,5	3	113	80058232



#### XRF Casting dishes suitable for Herzog (round version)

Diameter internal		Flange diameter (FD)	Height (H)	Base thickness (s)	Weight approx.	Heraeus drawing
top (OD)	bottom (UD)					
mm	mm	mm	mm	mm	g	
31	29	44	3,5	2,5	65	80308810
34	32	44	3,5	2,5	66	80623108
41	39	50	5,0	2,5	100	80623113



Europe

#### Heraeus Deutschland GmbH & Co. KG

Heraeus Performance Products  
 Heraeusstr. 12-14  
 63450 Hanau, Germany  
 Phone +49 6181.35-3740  
 Fax +49 6181.35-8620  
 heraeus-platinum-labware@heraeus.com  
 www.pt-labware.de

United States of America

#### Heraeus Precious Metals

**North America Conshohocken LLC**  
 Functional Materials  
 24 Union Hill Road  
 W. Conshohocken, PA 12-428-2736, USA  
 Phone +1 484 391 2527  
 Fax +1 610 825 7061  
 heraeus-platinum-labware@heraeus.com  
 www.pt-labware.com

