



XRF Platinum Program for Perl-X




XRF Crucibles suitable for Perl-X (cylindrical version with rim)

Diameter internal		Flange diameter (FD)	Height (H)	Base thickness (s)	Weight approx.	Heraeus drawing
top (OD)	bottom (UD)					
mm	mm	mm	mm	mm	g	
40	40	48	39	0,8	80	80623294
40	40	48	39	1,0	100	80623294




XRF Crucibles suitable for Perl-X (conical version)


Diameter internal		Flange diameter (FD)	Height (H)	Base thickness (s)	Weight approx.	Heraeus drawing
top (OD)	bottom (UD)					
mm	mm	mm	mm	mm	g	
40,4	38	50	39	0,8	78	80628592
40,4	38	50	39	1,1	98	80628590



XRF Lids suitable for Perl-X


Diameter	Total length (L)	Base thickness (s)	Weight approx.	Heraeus drawing
mm	mm	mm	g	
55	60	0,5	27	80068494





XRF Casting dishes suitable for Perl-X

Diameter internal top (OD)	Diameter internal bottom (UD)	Flange diameter (FD)	Height (H)	Base thickness (s)	Weight approx.	Heraeus drawing
mm	mm	mm	mm	mm	g	
32,4	30,0	65	3,6	0,8	45	80623282
32,0	30,5	55	2,5	1,0	50	80623292
32,0	30,5	65	3,5	1,0	53	80050321
35,5	33,5	65	3,4	0,8	48	80623284
37,0	34,5	55	3,5	1,0	46	80046700
40,5	39,5	60	3,5	1,0	65	80623287
40,5	39,5	65	3,5	0,8	50	80077497
40,5	39,5	65	3,5	1,0	60	80077497



Dispersion Hardened (DPH) Alloy

Through a special process, Heraeus has created a Dispersion Hardened (DPH) alloy.

The inclusion of finely dispersed zirconia offers an economical way to increase the mechanical strength of the product resulting in longer service times.

Engraving possibilities



For easier tracking of your platinum products in the lab Heraeus offers different engraving possibilities.

Europe

Heraeus Deutschland GmbH & Co. KG

Heraeus Performance Products
 Heraeusstr. 12-14
 63450 Hanau, Germany
 Phone +49 6181.35-3740
 Fax +49 6181.35-8620
 heraeus-platinum-labware@heraeus.com
 www.pt-labware.de

United States of America

Heraeus Precious Metals

North America Conshohocken LLC
 Functional Materials
 24 Union Hill Road
 W. Conshohocken, PA 12-428-2736, USA
 Phone +1 484 391 2527
 Fax +1 610 825 7061
 heraeus-platinum-labware@heraeus.com
 www.pt-labware.com

